

# Gluco Navii DMS user guide

Version 1.4.1.28 / 2014-05-21

## GlucoNavii DMS

The SD Diabetes Management Software is a personal computer (PC) software application that provides healthcare professionals to upload test data from meters, store the data in the Health Care System application and analyze the data by graphs and other software tools.

SD Diabetes Management Software provides you the following functions:

- Transmit blood glucose results from a meter to your PC.
- Transmit blood glucose results from your PC to E-mail and FTP server.
- Use electric data instead of a log book.
- Present the trends of test results.
- Share your data with your doctors.

## Minimum requirements

<b>CPU</b>	Intel Pentium (x86, x64) 800MHz
<b>RAM</b>	512MB
<b>HDD</b>	50MB
<b>I/O</b>	USB Controller
<b>Display</b>	1024 x 768
<b>OS</b>	Windows XP, VISTA, 7, 8

## Download at

<http://www.sdbiosensor.com/Program/5810>

## Installer

Setup.Product Name(Version Info).exe

## Components


- GlucoNavii DMS driver (default)
- SD USB Cable driver (optional )
- NFC R/W(Dragon) driver (optional)

## Install location

{System Directory}\GlucoNavii DMS\ (default)

## Starting GlucoNavii HCP

(Windows 7) start > All programs > GlucoNavii DMS > Click GlucoNavii DMS ICON.




<figure 1. – GlucoNavii DMS program>


## GlucNavii DMS interface

GlucNavii DMS interface is divided into 2 sections.

- User management section
- Data management section




<Figure 2 – User management section screen>



<Figure 3 – Data management section screen>


## GlucoNavii DMS Setting

Click “About” button in the navigation part – you are able to check version of DMS, copyright and homepage link.




<Figure 4 – About>

Click “Preferences” button in the navigation part – you are able to modify unit, glucose target range, language and passcode.



<Figure 5 – Preferences>


### Preference- glucose

A screenshot of the 'Preferences' dialog box. The 'Glucose' tab is selected, and the 'Localization' tab is also visible. The 'Unit' is set to 'mmol/L'. The 'Hypo' value is 70.0 and the 'Hyper' value is 200.0. There are 'Ok' and 'Cancel' buttons at the bottom. The dialog box has a close button 'X' in the top right corner.

Field	Value
Unit:	mmol/L
Hypo:	70.0
Hyper:	200.0

<Figure 6 – Glucose>

## Preferences – localization




The screenshot shows a 'Preferences' dialog box with a close button (X) in the top right corner. It has four tabs: 'Glucose', 'Localization', 'Account', and 'Graph'. The 'Localization' tab is selected. Inside the dialog, there are two settings: 'Language' is set to 'English' with a small US flag icon to the left of the text; 'Date format' is set to 'yyyy/mm/dd'. At the bottom of the dialog are two buttons: 'Ok' and 'Cancel'.

<Figure 7 – Localization>

## Preferences – account

GlucoNavii DMS provides one account. When you check “application lockout” and set up the password, you are not able to start program without entering the password.




The screenshot shows the same 'Preferences' dialog box, but with the 'Account' tab selected. The 'Option' section has a checkbox labeled 'Application lockout' which is currently unchecked. Below this are two text input fields: 'Password:' and 'Re-enter password:'. At the bottom of the dialog are two buttons: 'Ok' and 'Cancel'.


<Figure 8 – Account>

## Preferences – graph

- Y-Axis – when “Auto range” is activated, the scale of the vertical values are changed automatically in the chart depends on the glucose level. For example, if you’re maximum glucose level is 300, the maximum value of Y-axis will be 300.
- Layout – there are two types of layout you can choose including “single view” and “tap view”. The default layout is single page view.



<Figure 9 – Graph>



<Figure 10 – Tap view>


## Connecting a meter to your PC

Connecting method – you are able to transfer the data from the meter to the PC through either NFC RW or USB cable. If NFC RW and USE cable are connected simultaneously to the PC, only one device which is connected first will transfer the data to the PC.

- NFC RW↔PC connecting: tap the back of the meter on the NFC RW for 1~ 2 sec.

Make sure that you wait for 1~2 sec in order to transfer the data perfectly


- USB Cable↔PC connecting: you must connect the cable between PC and meter when the meter is off. After connecting, you can turn on the meter to transfer the data



<Figure 11– Competing data transfer>


## Creating new user or adding to existing user

After completing data transfer, the data is automatically sorted by serial number. If there is no previous serial number, you are able to either create new user or add it to existing user.



<Figure 12– Creating new user>

Click “Create New User” button – you are able to create new user and enter the user’s information



The image shows two side-by-side screenshots of the 'User Information' dialog box. The left screenshot shows the 'Create New User' path, and the right screenshot shows the 'Add to existing user' path. Both screenshots show the same form fields, but the right one has pre-filled values for 'First name' (Users), 'Last name' (SD), 'Birthdate' (1985-07-07), 'Height' (165 cm), 'Weight' (60 kg), and 'Diabetic since' (2012-05-12).

**User Information** ×

First name  
Last name  
e-mail  
Gender  
Female  
Birthdate  
2013  
Height  
cm  
Weight  
kg  
Diabetic type  
Insulin dependent  
Diabetic since  
2013  
Target blood glucose  
70 - 203 mg/dL  
Device list  
No Serial no  
New M07A01AA9007  
Ok Cancel


**User Information** ×

First name  
Users  
Last name  
SD  
e-mail  
Gender  
Female  
Birthdate  
1985 7 7  
Height  
165 cm  
Weight  
60 kg  
Diabetic type  
Insulin dependent  
Diabetic since  
2012 5 12  
Target blood glucose  
70 - 203 mg/dL  
Device list  
No Serial no  
New M07A01AA9007  
Ok Cancel

<Figure 13 – User Information>




After completing user registration, transferred data will be appeared (figure 14).




<Figure 14 – Data management screen>

Click “add to existing user” button – you are also able to add the data to the existing user



<Figure 15 – Add to existing user>



<Figure 16 – Add to existing user dialog>

Select the user you want to add, and then click "OK" button.

## Log book

Column group panel			
Date	Time	Event	
10/14/2013: 1			
10/14/2013	1:20 PM		148
9/25/2013: 3			
9/25/2013	7:15 PM		Lo
9/25/2013	7:14 PM		41
9/25/2013	7:13 PM		55
7/25/2013: 15			
7/25/2013	5:56 PM		108
7/25/2013	5:56 PM		103
7/25/2013	5:53 PM		101
7/25/2013	5:52 PM		99
7/25/2013	5:51 PM		97
7/25/2013	5:49 PM		95
7/25/2013	5:48 PM		100
7/25/2013	5:38 PM		285
7/25/2013	5:37 PM		308
7/25/2013	5:35 PM		293
7/25/2013	5:33 PM		294
7/25/2013	5:32 PM		283
7/25/2013	5:30 PM		287
7/25/2013	5:29 PM		279
7/25/2013	5:27 PM		279

<Figure 17 – Logbook grid part >

The data is grouped by date and displayed in descending order.

### <Removing column group>


Drag and drop the column which you want to remove from “column group panel” to “column heater”



<Figure 18 – Removing column group>


### <Grouping column>


Drag and drop the column which you want to group by from “column heater” to “column group panel”




<figure 19 – Grouping column>



### <Sorting column>

Click  icon after placing the mouse pointer to the column header- you are able to sort the column in descending/ ascending order.




### <Filtering column>

Click  icon after placing the mouse pointer to the column header- you are able to select glucose or event filter options (figure 18&19)




<Figure 10– Glucose filter option >



<Figure 21– Event filter option >

## Graphs



<Figure 22 – Graph >

This graph shows glucose results for a selected date range. This graph also displays glucose target range, hyperglycemia and hypoglycemia in order to help you easily track changes in glucose results. .

### <Change glucose target range>

Click “target range” button – you are able to change the glucose target range

If the user is previously registered, <figure 23> will be appeared, if not, <figure 24> will be appeared

**User Information** ×

First name

Last name

e-mail

Gender

Birthdate

Height   Weight

Diabetic type

Diabetic since

Target blood glucose  
 -  mg/dL

Device list

No	Serial no
1	M07A01AA9007

<Figure 23 – Changing glucose target range – registered user>

**Preferences** ×

Glucose Localization Account Graph

Unit:


Hypo:

Hyper:


< Figure 24 – Changing glucose target range – not registered user>

**<Selecting the date interval>**

Click “date interval” button <figure 25>- you are able to select the date interval 7/14/30/ 90 days or customize the days you want.



<Figure 25 – Click data interval button>



<Figure 26 – Selecting date interval >



## Statistics

	Avg	Min	Max	SD	Count of measurement			
					Hypo	Target	Hyper	Sum
					Analysis period: Total: 2013/06/01 ~2013/06/11 <span style="float: right;">300</span>			
<b>Pre-meal</b>	13.9	5.8	22.1	4.9	0	48	67	115
<b>Post-meal</b>	14.0	5.9	22.2	4.8	0	45	68	113
<b>Other</b>	9.1	5.6	21.3	4.8	0	33	39	72
<b>Total</b>	13.9	5.6	22.2	4.8	0	126	174	300
<b>Mon</b>	13.7	8.3	19.2	4.9	0	24	24	48
<b>Tue</b>	12.9	9.6	19.8	4.3	0	24	10	34
<b>Wed</b>	11.5	10.9	12.2	0.4	0	24	0	24
<b>Thu</b>	12.9	12.3	13.5	0.4	0	6	18	24
<b>Fri</b>	14.2	13.6	14.9	0.4	0	0	24	24
<b>Sat</b>	14.3	5.6	21.1	5.7	0	24	54	78
<b>Sun</b>	15.1	6.9	22.2	5.9	0	24	44	68

<Figure 27 – Statistics>

The glucose results with control solution are excluded from statistic figures. The number of values with control solution are displayed in the top right corner of a screen.

	Avg	Min	Max	SD	Count of measurement			
					Hypo	Target	Hyper	Sum
					Analysis period: Total: 06/01/2013 ~06/11/2013 <span style="float: right;">300 297 3</span>			
<b>Pre-meal</b>	251	104	398	87.9	0	37	78	115
<b>Post-meal</b>	252	107	400	87.0	0	36	77	113
<b>Other</b>	150	101	383	84.1	0	27	42	69
<b>Total</b>	249	101	400	86.5	0	100	197	297
<b>Mon</b>	247	149	346	88.2	0	24	24	48
<b>Tue</b>	233	173	356	77.5	0	24	10	34
<b>Wed</b>	208	197	220	7.1	0	4	20	24
<b>Thu</b>	232	221	244	7.1	0	0	24	24
<b>Fri</b>	256	245	268	7.1	0	0	24	24
<b>Sat</b>	253	101	380	103.5	0	24	51	75
<b>Sun</b>	272	125	400	106.5	0	24	44	68

<Figure 28 – Number of test result with control solution>


### <Changing analysis period>


Statistic table shows each categories based on data in logbook. Thus you are able to see the statistic for certain period by editing the analysis period

Analysis period:	Total: 3/10/2013 ~10/14/2013				
	Total: 3/10/2013 ~10/14/2013				
	Graph: 10/14/2013 ~10/14/2013				
	<b>SD</b>	<b>Count of measurement</b>			
		<b>Hypo</b>	<b>Target</b>	<b>Hyper</b>	<b>Sum</b>

<Figure 29 – Changing analysis period >

## Printing

Click  button in user information part – you are able to print out the report, save the data and attach the data to an email.



The screenshot shows the "Export" dialog box in the SD Biosensor reporting software. The dialog is titled "Report" and contains the following sections:

- Report:** Includes a "Print" button, a "Copy" dropdown set to "1", "Printer properties", a "Save as PDF" dropdown, and a "Quality" slider set to 90.
- Report Settings:** Includes an "Edit Remark" button and "Header / Footer Settings".
- Preview Settings:** Includes a zoom slider set to "x 1.00", radio buttons for "Page width" (selected) and "Whole page", and a "Close" button.

The main preview area displays the report content for "User SD":

- SD BIOSENSOR, INC. logo and name.
- Birthdate: 7/7/2013 ( Female )
- Date of service: 1/17/2014
- Glucose graph ( 10/14/2013 ~ 10/14/2013 )
- Target range: Hypo: 70, Hyper: 300
- Analysis matrix ( 10/14/2013 ~ 10/14/2013 )

	Avg	Min	Max	SD	Measurement count			Sum
					Hypo	Target	Hyper	
Pre-meal	0	0	0	0,0	0	0	0	0
Post-meal	0	0	0	0,0	0	0	0	0
Other	148	148	148	0,0	0	1	0	1
Total	148	148	148	0,0	0	1	0	1
Mon	148	148	148	0,0	0	1	0	1
Tue	0	0	0	0,0	0	0	0	0
Wed	0	0	0	0,0	0	0	0	0
Thu	0	0	0	0,0	0	0	0	0
Fri	0	0	0	0,0	0	0	0	0
Sat	0	0	0	0,0	0	0	0	0
Sun	0	0	0	0,0	0	0	0	0


<Figure 30 – Printing dialog>

### <Printing>

Click "print" button to print out the report

### <Changing printer>

Click ▼ (dropdown menu) next to the print button in order to change the printer




<Figure 31 – Selecting printer>

### <Saving / attaching file>

You are able to save the report as PDF/JPG file or attach the report to send email through Microsoft outlook.


Click ▼ (dropdown menu) next to the save as button




<Figure 32 – Saving and sending the report >


<Editing note>

Note section represents the red box in the report as below




<Figure 33 – Note section>

Click  button or note section directly in the report preview in order to edit the note.



<Figure 34 – Editing note >




<Figure 35 – Click edit remark >

<Editing header /footer>

Click “header/footer settings” to edit header and footer

The screenshot shows the 'Export' window. On the left, there are three sections: 'Report', 'Report Settings', and 'Preview Settings'. In the 'Report Settings' section, the 'Header / Footer Settings' button is highlighted with a red box. The main preview area shows a report for 'Gemma Leem' with a birthdate of 6/5/1979 and a date of service of 5/13/2014. It includes a 'Glucose graph ( 2/3/2014 ~ 2/9/2014 )' and an 'Analysis matrix ( 2/3/2014 ~ 2/9/2014 )' table.

	Avg	Min	Max	SD	Measurement count			
					Hypo	Target	Hyper	Sum
Pre-meal	218	50	245	63.2	1	0	8	9
Post-meal	73	65	90	8.5	2	5	0	7
Other	67	60	247	72.6	10	21	9	40
Total	120	50	247	70.1	13	26	17	56
Mon	108	50	234	63.9	1	6	1	8
Tue	119	80	231	65.7	2	5	1	8
Wed	131	61	245	68.6	3	2	3	8
Thu	131	61	243	68.0	3	2	3	8
Fri	136	66	246	69.7	2	3	3	8
Sat	136	65	245	67.8	1	4	3	8
Sun	136	66	247	67.1	1	4	3	8




<Figure 36 – Editing header /footer>

You are able to choose the header or footer by selecting tab on the top of the screen

#### <Editing background>

You are able to select either photo (BMP, JPEG, PNG) or single color for background of header and footer




<Figure 37 – Selected monochrome for header >

#### <Editing sentence in header/footer>

After checking  Display ”, you are able to edit font size, color, and style in header/ footer.







<Figure 38 – Example of sentence in header>

<Reset header/footer>

Click  button in order to reset the header/footer

## Exporting

Click  or  button in user information part – you are able to export the data into excel /CSV format and attach the data to an email.



<Figure 39 – Exporting icons (excel and CSV)>

<Exporting glucose data into an excel/CSV format>

**Export glucose data** [X]

**Output File**

Save to File  
C:\Users\Wjhjang\Documents\GlucoNaviiDMS (5).xls [...]

Send an Email

**Data Selection**

Graph Period (34, 06/10/2013 ~ 06/11/2013)

Complete Period (300, 06/01/2013 ~ 06/11/2013)

Custom Period (300, 06/01/2013 ~ 06/11/2013)  
06/01/2013 ~ 06/11/2013

**Output Format**


Excel

CSV

OK Cancel


<Figure 40 – Export glucose data>

After choosing the location where you want to save the file and selecting the file format, click OK button in order to export the glucose data into an excel/ CSV format




< Figure 41 – Operation complete >

<Attaching the file and send an email>



< Figure 42 – Send an Email >

## User management


Click  icon next to the “Export” button – you are able to move to user management screen







<Figure 43 – Icon for user who is not registered >




<Figure 44 – Icon for user who is registered>




<Figure 45 – User management screen >


Several function icons are displayed next to the each user name (data  / export  / edit  / delete  )

### <Checking user information>

Click  icon, and then user information will be appeared the right side of the screen

About		Preferences	
		Sort by ▼	
Data Count	Last Update		
227	10/16/2013		

<Figure 46 – User information icon >




The screenshot shows the GlucoNavii DMS interface. At the top, there are 'About' and 'Preferences' icons. Below them is a search bar and a 'Sort by' dropdown menu. The main content area displays a table with columns for 'Name', 'Data Count', and 'Last Update'. The table contains one entry: 'UserSD' with a data count of 227 and a last update date of 10/16/2013. To the right of the table, there is a 'UserSD' information panel with fields for 'e-mail', 'Gender' (Female), 'Birthdate' (7/7/2013), 'Height' (165.0 cm), 'Weight' (60.0 kg), 'Diabetic type' (Insulin-Dependent), 'Diabetic since' (5/12/2013), and 'Target glucose' (70 - 203 mg/dL). Below this is a 'Device list' section with the entry 'M07A01AA9007'. The footer of the interface includes 'SD Biosensor, Inc.' and icons for 'Home' and 'YouTube'.

<Figure 47 – User information>

### <Sorting user column >

Click “sort by” button – you are able to sort the user column by last update/name/ data count

Sort by ▼		
a Count	<u>L</u> ast Update	
2	<u>N</u> ame	
	✓ <u>D</u> ata Count	

<Figure 48 – Sort by>