OBSAH

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FDIR-V12

PN:6004112F001-03A

The production date is shown on the outer package.

For American please refer to °F,
For European please refer to °C

Before using your device, please read these instructions and warnings completely.

CE 1639

Description Version: 022025

EN PRODUCT INSTRUCTIONS

CONTENT

1 Product Description	31
2 Product structure and components	31
3 Product specifications	31
4 Intended use	32
5 Contraindications	32
6 Notes	32
7 Installation and use	34
8 Body Temperature	36
9 Troubleshooting instructions	37
10 Maintenance and maintenance	37
11 Packing List	38
12 Final disposal	38
13 Explanation of standardized symbol	38
14 Electromagnetic compatibility information	39
15 Warranty	42


1. Product Description

Thank you very much for choosing our products! This high-tech infrared thermometer is used to measure human body temperature by absorbing infrared energy from forehead of human. It can help you learn about the health of you and your family quickly at any time and anywhere.
Product name: Infrared thermometer
Product model: FDIR-V12


2. Product structure and components

The thermometer is mainly composed of plastic shell, infrared temperature sensor, PCBA, buzzer, LCD display screen and battery.

Main appearance structure diagram description:

Display Description:

3. Product specifications

Product name	Infrared thermometer
Product model	FDIR-V12
Measuring mode	Forehead/Object Surface
Reference site	Oral cavity
Power supply	DC 3V, 2x 1.5V AAA batteries
Measurement range	Body temperature: 32.0°C~42.9°C (89.6°F~109.2°F) Object temperature: 0°C~60.0°C (32.0°F~140.0°F)
Measurement accuracy	Body temperature: ±0.2°C between 35.0°C~42.0°C (±0.4°F between 95.0°F~107.6°F) ±0.3°C (or ±0.5°F); outside this ranges Object temperature: ±1°C (1.8°F)
Digital increment of displayed temperature	0.1°C (0.1°F)
Clinical repeatability	Within ±0.3°C (0.5°F)
Measurement time	1s
High body temperature hint	≥37.6°C (100.4°C)

Operation environment	Temperature: 10°C~40°C (50°F~104°F) Relative humidity: 15%~95% RH, No condensation Atmospheric pressure: 70kPa~106kPa
Transportation and storage environment	Temperature: 25°C~55°C (-13.0°F~131.0°F) Relative humidity: 15%~95% RH, No condensation Atmospheric pressure: 70kPa~106kPa
Temperature unit	°C/°F
Measurement memory	19 the memories
Product size	42 mm × 52 mm × 144 mm
Net weight	70 g (Battery free)
Period of validity	5 years
Software version	V1.0
Grade of waterproof	IP22
Electric shock	Internally powered ME equipment
Operational mode	Continuous operation
Applied part	Type BF applied part, including the whole unit.

*Statement: If stored or used outside the temperature, humidity or atmospheric pressure range specified by the manufacturer, the product may not achieve the claimed performance

4. Intended use

1. The product displays the body temperature of the forehead through thermal radiation.
2. Expected patients: and can be used for body temperature measurement in infants, children and adults.
3. Intended user: Doctor, nurse or regular user. (We recommend adults to operate the thermometer instead of babies and children.)
4. Expected Operating Environment: Hospitals, clinics and home Settings

5. Contraindications

Do not use local lesions such as inflammation, trauma and postoperative.

6. Notes

- 6.1 Measurement
 1. This product is only for monitoring and self-testing body temperature, can not be used for disease diagnosis, patients only through the results of self-judgment and treatment is very dangerous, so please follow the guidance of the doctor.
 2. To ensure the accuracy of the measurement, please set the thermometer in the normal working environment for more than 30 minutes, so that it reaches the equilibrium state before measuring.
 3. Before measurement, ensure that the probe lens surface is not covered with small foreign bodies (such as dander, droplets, dust, etc.), so as not to affect the measurement results.

4. Please do not measure under the condition that the ambient temperature changes rapidly, such as the air outlet of the air conditioning or heater, which may affect the measurement results.
5. Please do not measure under strong electromagnetic interference environment (working microwave oven, induction cooker, mobile phone nearby, etc.), this may cause measurement is not correct or can not be measured.
6. Please ensure that there is no sweat, cosmetics, oil stains on the forehead of the test subject, so as not to affect the accuracy.
7. When the subject enters the warm room from the cold outside, the forehead temperature will be temporarily lower than the real temperature. Please wait for 30 minutes before measuring.
8. When the subject's forehead is covered with long hair, the temperature of the forehead will be temporarily higher than the real temperature. Please wait for 30 minutes before measuring.
9. Participants were measured after taking medications that change their body temperature (e.g. aspirin, acetaminophen, ibuprofen), which may cause measurement bias. 10. The subjects may cause temperature fluctuations after intense exercise, crying, eating, etc. It is recommended to keep quiet for 30 minutes before measurement.
11. When measuring, please keep the subject quiet and steady, according to the measurement method shown in the specification, so as not to cause the measurement deviation caused by improper operation.
12. When measuring, please try to aim at the center of the forehead and keep the specified measuring distance with the probe to avoid measurement deviation.
13. Continuous measurement, the measurement results may have a small error, which is a normal phenomenon, it is recommended that the maximum continuous measurement in a unit of time 3 times.
14. Different types and brands of temperature measuring products do not have comparative value. It is recommended to use the same temperature measuring product to observe temperature changes in order to better grasp the temperature changes.

6.2 About Products

1. The infrared probe is the most easily damaged part of this product. Please avoid touching the probe with your finger or blowing the probe with your mouth, which may cause measurement deviation.
2. Before use, confirm whether the probe lens is dirty, if dirty should follow the following steps to clean, otherwise it will cause measurement deviation:- Gently wipe with cotton swabs or soft cloth dipped in more than 75% concentration of medical alcohol;
3. Please do not force collision, fall, stamped or shake this product, may cause product failure or cause measurement deviation.
4. Do not disassemble, repair or modify the product.
5. Do not use in electromagnetic interference environment, which may affect the measurement results.

6. Please keep this product out of reach of children to prevent children from swallowing batteries or small parts. If swallowed, please contact the doctor immediately.
7. Please keep the product in accordance with the product transportation / storage conditions to avoid use problems or measurement deviations.
8. Do not put in direct sunlight, high temperature and humidity, dust, near the fire, vulnerable to vibration shock place custody.
9. The thermometer is not used for a long time, the battery should be removed to prevent battery leakage.
10. Do not place this product in a place with electric shocks.
11. Do not throw the thermometer and its accessory batteries into the fire, which may cause an explosion.
12. Please dispose of waste or residue at the end of the product's service life in accordance with local laws and regulations.
13. The normal use of this product is the human body temperature measurement mode. If you need to enter the correction mode must contact the company's after-sales service agencies.
14. In case of any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the member state.


7. Installation and use

7.1 Installation of batteries

Put two batteries in the packaging box into the battery house on the back of the device. At the point the product will start self-inspection, and then will be ready to start measurement. (If the battery power is low when starting up, please replace the battery to ensure adequate power supply)

7.2 The forehead temperature measurement


1. Aim the thermometer probe at the center of the forehead and keep a distance of 1 cm~6 cm.
2. Press "measuring key", the thermometer will start up and automatically measure the temperature.
3. About 1s later, the thermometer displays the measurement result with a prompt sound. Note: Forehead mode will be entered by default every time you start the machine, and "°C" forehead icon will be displayed in forehead temperature mode.

7.3 Object mode temperature measurement

1. Tap "Mode key" to switch to the object temperature mode. "°F" Object icon will be displayed in object Temperature mode.

2. Align the thermometer with the object and press the "measuring key" button to measure the temperature. After about 1s, the thermometer displays the measurement results with a prompt sound.

7.4 Display and Sound Tips

1. After the measurement is completed, the following tips may appear in the thermometer:

Body Mode Screen display	Sound Tips	Back light color	Representation of information
Temperature: 32.0°C~37.5°C (89.6°F~99.65°C)	"beep~"	green	Target temperature values for this measurement
Temperature: 37.6°C~42.9°C (99.7°F~109.2°C)	"beep beep beep~" "beep beep beep~" "beep beep beep~"	red	Target temperature values for this measurement
Err	"beep~" "beep~" "beep~"	green	The ambient temperature is not within the allowable operating range of the thermometer
Hi	"beep~" "beep~" "beep~"	green	The measurement temperature is higher than the maximum value of the thermometer range
Lo	"beep~" "beep~" "beep~"	green	The measurement is below the minimum value of the thermometer range

2. Thermometer presses the "measuring key" again within 30 s after the measurement is completed, and the next measurement can be carried out.

7.5 Automatic shutdown function

If the thermometer is not operated within 30 s, it will automatically enter the shutdown state.

7.6 Memory query function


When the thermometer is on, tap "Memory/Voice key" to enter the memory query interface. The screen displays the last effective temperature measured. Tap "Memory/Voice key" again, you can query the previous effective measured temperature value, and so on, can be recycled to query up to 19 memory value.

7.7 Changing the Measurement Unit

1. In the boot state, hold down the "Mode key" and enter the temperature unit setting interface after about 10 s.
2. keep long press "Mode key", or units automatically switch every 2 s, when the desired unit is displayed, release

the hand, that is, select the current display unit, the thermometer immediately into the state to be measured.

7.8 Mute setting function


- After power on, press [Memory/Voice Button] for about 3 seconds to turn on or off the mute mode.
- When the mute function is enabled, an "M" icon is displayed on the screen. When the mute function is off, an "M" icon is displayed on the screen.

7.9 Low Voltage Prompt Function

- The thermometer automatically detects the battery voltage every time it is turned on. When the battery voltage is insufficient but can continue to use, the "Low Voltage" icon appears on the display screen, but does not affect the operation.
- When the battery voltage is too low to continue to use, the "LO" icon is displayed and the "Low Voltage" icon blinks, about 5 s after forced shutdown, thermometer needs to replace the new battery to continue to use.

8. Body Temperature

The normal body temperature of people is a range value, the normal body temperature of different people will be different, and the personal body temperature will also change at different times. The following are the reference values for the definition of fever in various parts of the human body:

- *Data sources: ISO 80601-2-56:2017/AMD 1:2018, for information only.
- *This infrared thermometer has the same forehead temperature measurement pattern as the sublim temperature. It is dangerous to make self-diagnosis and treatment according to the measurement results. Please consult a professional doctor with the measurement results.

9. Troubleshooting instructions

Fault Description	Possible causes	Solution
The screen does not react when the power is on or equipped with batteries.	Battery runs out Battery installation error	Please replace the new battery and re-open Please install the battery correctly according to the positive and negative pole mark in the battery bin.
	After reloading the battery, the screen remains blank	Please contact after sale
	The ambient temperature is not within the allowable operating range of the thermometer	Please place the thermometer in a room with normal working environment temperature and keep it for 30 minutes before measuring.
	Hardware failure	Please contact after sale
	Too low battery voltage, Causing failure to boot.	Please replace the new battery and restart
	Battery power is running out	Please replace the new battery as soon as possible
	The measurement temperature is higher than the maximum value of the thermometer range	Please check the measurement target is correct and re-measure again.
	The measurement temperature is below the minimum value of the thermometer range	Please check the measurement target is correct and re-measure again.

10. Maintenance and maintenance

- Please use a soft, dry cloth to gently wipe the dirt on the mainframe. If the surface is more serious dirty, you can dip some alcohol to wipe. Cleaning is recommended every 1 month.
Caution: Do not wash with water or use cleaners, diluents or volatile oils containing abrasives.

- After the thermometer probe is dirty, please gently wipe it with cotton swabs or soft cloth dipped in more than 75% concentration of medical alcohol. If it is still not clean, please contact customer service. Note: do not use paper towels directly, may scratch thermometer probe lens.
- If the thermometer probe or key, display and other parts damaged, please contact after-sale processing.
- Do not put in direct sunlight, high temperature and humidity, dust, near the fire, vulnerable to vibration shock place custody.
- When the thermometer is not used for a long time, the battery should be removed to prevent battery leakage.
- When the "batteries icon" appears on the thermometer display screen, it indicates that the battery power is about to run out. Please replace the battery as soon as possible.
- We do not authorize any organization or individual to carry out maintenance, please do not disassemble and adjust the thermometer.
- The thermometer has been tested and calibrated at the time of manufacture. If you have any questions about the function or measurement accuracy of the product, please contact the dealer or manufacturer.

11. Packing List

Use only original accessories to check whether the received accessories are complete.

Quantity	Components
1 pcs	FDIR-V12 host
2 pcs	AAA batteries
1 pcs	User Manual

*The technical data should be combined with the user manual.

*In addition to specifying the use of AAA batteries, replacing existing components with non-manufacturer supplied components may cause measurement errors and reduce safety.

12. Final disposal

At the end of the product life cycle, do not throw this product into the normal household garbage, but bring it to a collection point for the recycling of electronic equipment.

Waste Electrical and Electronic Equipment can have potentially harmful effects on the environment. Incorrect disposal can cause harmful toxins to build up in the air, water and soil and can be harmful to human health.

13. Explanation of standardized symbol

CE 1639	Complies with the European Medical Device Regulation (EU) 2017/745, Notified Body is SGS Belgium NV.
EU REP	Authorized representative in the European Community.
	Attention: see Instructions for use!

	Caution! Consult accompanying documents.
	Type BF applied parts
LOT	Batch code
	Manufacturer information:
IP 22	IP code of the device: this device's grade of against ingress of solid foreign objects.
	Disposal in accordance with Directive 2002/96/EC (WEEE).
MD	Medical Device
UDI	Unique Device Identifier
	Temperature limit
	Humidity limitation
	Atmospheric pressure limitation
	Date of manufacture

14. Electromagnetic compatibility information

WARNING:
Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Not use of accessories, transducers and cables other than those specified or provided by the manufacturer of this the FDIR-V12 could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the FDIR-V12 including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Guidance and manufacturer's declaration – electromagnetic emission – for all EQUIPMENT AND SYSTEMS

Guidance and manufacturer's declaration – electromagnetic emission		
The FDIR-V12 is intended for use in the electromagnetic environment specified below. The customer or the user of FDIR-V12 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR11	Group 1	The FDIR-V12 uses RF energy only for its internal function. There for, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

RF emissions CISPR11	Class B	The FDIR-V12 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	N/A	

Guidance and manufacturer's declaration – electromagnetic immunity – for all EQUIPMENT AND SYSTEMS

Guidance and manufacturer's declaration – electromagnetic immunity			
The FDIR-V12 is intended for use in the electromagnetic environment specified below. The customer or the user of the FDIR-V12 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be covered with synthetic material, the relative humidity should be at least 30 %. If ESD interference with the operation of equipment, counter measurements such as wrist strap, grounding shall be considered.
Electrostatic transient/burst IEC 61000-4-4	± 2 kV for power supply lines	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT (100 % dip in UT) for 0.5 cycle 0 % UT (100 % dip in UT) for 1 cycles 70 % dip in UT for 25/30 cycles 0 % UT (100 % dip in UT) for 250/300 cycles	N/A	Mains power quality should be that of a typical commercial or hospital environment. If the user of the FDIR-V12 requires continued operation during power mains interruptions, it is recommended that the FDIR-V12 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity – for EQUIPMENT AND SYSTEM that are not LIFE-SUPPORTING

The FDIR-V12 is intended for use in the electromagnetic environment specified below. The customer or the user of the FDIR-V12 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the FDIR-V12, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: PORT IMMUNITY TO RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)
Radiated RF IEC 61000-4-3	6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 10 V/m 80 MHz to 2.7 GHz 385 MHz-5785 MHz Test specifications for ENCLOSURE PORT IMMUNITY TO RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)	10 V/m 80 MHz to 2.7 GHz 385 MHz-5785 MHz Test specifications for ENCLOSURE PORT IMMUNITY TO RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)	$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.

- The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz, 13.553 MHz to 13.567 MHz, 26.957 MHz to 27.283 MHz, and 40.66 MHz to 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.
- Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed

RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the FDIR-V12 is used exceeds the performance is observed, additional measures may be necessary, such as reorienting or relocating the FDIR-V12.
c Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM - for EQUIPMENT AND SYSTEMS

Recommended separation distances between portable and mobile RF communications equipment and the FDIR-V12				
The FDIR-V12 is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the FDIR-V12 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the FDIR-V12 as recommended below, according to the maximum output power of the communications equipment.				
Rated-maximum output of transmitters (W)	Separation distance according to frequency of transmitter (m)			
	150 kHz to 80 MHz outside ISM and amateur radio bands $d = \left[\frac{1}{f} \right] \sqrt{P}$	150 kHz to 80 MHz in ISM and amateur radio bands $d = \left[\frac{1}{f} \right] \sqrt{P}$	800 MHz to 2.7 GHz $d = \left[\frac{1}{f} \right] \sqrt{P}$	800 MHz to 2.7 GHz $d = \left[\frac{1}{f} \right] \sqrt{P}$
0.01	0.12	0.20	0.035	0.07
0.1	0.38	0.63	0.11	0.22
1	1.2	2.00	0.35	0.70
10	3.8	6.32	1.10	2.21
100	12	20.00	3.5	7.0

For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

15. Warranty

We provide one year warranty starting from the date of purchase. Please refer to the followings situations that are excluded from the free repair services within the warranty period.
1. All damages caused by disassembly and repair of the device by yourselves.
2. All damages caused by dropping the device during usage, or transport.
3. All damages caused by improper usage of the device and not following the instructions on the user manual.

Please contact after-sales service and support and enclose your product purchase receipt while claiming for warranty services.

Location of purchase:

Contact number:

Date of purchase:

Famidoc Technology Co., Ltd.
No. 212 Yilong Road, Chang'an Town, Dongguan, Guangdong Province, 523853, P.R. China

EU REP
Name: Shanghai International Holding Corp. GmbH (Europe)
Add: Eilestrasse 80, 20537 Hamburg, Germany