



HEM-907 Extended Tube
Using the longer replacement tube (1.3 m length) allows your patient to sit back and relax.



HEM-907 Pole Mount
Bring the OMRON HEM-907 closer to the patient's bed, using this pole mount. Will also allow fixing to a bedrail. Including cuff storage basket.



HEM-907 Wall Holder
Fix the OMRON HEM-907 close to the patient using this wall holder with storage basket for cuff(s).




HEM-907 Cuffs To assure the highest versatility in using the HEM-907 monitor, OMRON offers these durable and wipeable cuffs in three sizes. It is essential to always use the appropriate size on your patients in order to achieve the highest accuracy.

- » OMRON CWS • 17–22 cm
- » OMRON CWM • 22–32 cm
- » OMRON CWL • 32–42 cm

COMPARISON TABLE

Blood Pressure Monitors Fully Automatic Upper Arm

| Model | ABPM 24/7 | HBP-1300 | i-Q132 | i-Q142 |
|-------------------------------------|----------------|--------------|-------------------------------------|------------------------------------|
| Model Number | BP5 | HBP-1300-E | HEM-1010-E | HEM-1040-E |
| Intellisense™ Technology | | | ✓ | ✓ |
| AC-Adapter | | ✓ | E1000 (incl. 3094100-9) | E1600 (incl. 3094298-6) |
| PC Connection with Software & Cable | ✓ | | | ✓ |
| Supplied Cuff | Medium/Large | Medium/Large | Small/ Medium arm in style 17–32 cm | Medium/Large arm in style 22–42 cm |
| Optional Large Cuff Available | ✓ | ✓ | | Built-in |
| Optional Small Cuff Available | ✓ | ✓ | Built-in | |
| Memory Positions | 400 | | 84 p. user 2 users | 84 p. user 2 users |
| Three-fold Displays (Sys/Dia/Pulse) | ✓ | ✓ | ✓ | ✓ |
| Date/Time Indication | ✓ | | ✓ | ✓ |
| Large Display and Operation Buttons | | ✓ | ✓ | ✓ |
| Dual Check System | | | | |
| Cuff Wrapping Guide | | | | |
| Hypertension Indicator | | | Symbol | Symbol |
| Irregular Pulse Detection | | | ✓ | ✓ |
| Averaging of 3 Values | Using Software | | ✓ | ✓ |
| Clinical Validation | IP | IP | IP | IP |



BLOOD PRESSURE MONITORS PROFESSIONAL

OMRON provides the gold standard for blood pressure monitoring in clinics and hospitals around the world. We've worked with doctors and clinicians to understand their needs and designed and modified our monitors accordingly. Exact results are guaranteed every time, that's why our monitors are used in 50% of clinical trials and national health surveys worldwide.

Our three key priorities in the professional field are accuracy, durability and functionality. These machines can be used hundreds of times a day, every day. They are sturdy to cope with the workload, while guaranteeing long-term reliability. And they do it quickly and simply.



Intelligent sensing helps people monitor and control their blood pressure so they can lead a healthy life without having to worry about inaccurate results or uncomfortable measurements.



BPS

ABPM 24/7

The Omron M24/7 ambulatory blood pressure monitor is a silent and compact device which can be easily taken anywhere with a patient ensuring maximum compliance. Using just 2AA batteries the M24/7 operates up to 48 hours and can store more than 400 readings in its solid-state memory for an unlimited length of time. Additional measurements can be triggered and specific events can be marked manually. The M24/7 can be connected to practically any standard PC with its optical USB interface cable. The device comes with two latex-free cuffs in medium and large wipeable cuff sleeves. The M24/7 has been clinically validated to the BHS and AAMI accuracy criteria.



HBP-1300-E

HBP-1300

The HBP-1300 Professional Blood Pressure Monitor is easy to use with a durable design that includes a built-in handle making the monitor easily portable. It is designed for use in professional settings and is clinically proven to produce fast and reliable results.

The monitor is compatible with a series of wipe-clean GS cuffs ranging from SS (12-18cm) to XL (42-50cm). It comes with an AC adapter as well as a rechargeable battery pack.

The device is fitted with a shock-proof bumper to protect it from accidental damage and comes with both a manual and an auscultatory mode. It is possible to manually adjust the inflation pressure for improved patient comfort and the Accuracy Check function includes regular pulse detection which displays an icon on the large LCD screen if the patient has an irregular heartbeat which would normally result in an inaccurate reading. In addition the Zero Indicator allows you to be confident that the monitor has been reset to zero allowing for a more accurate reading.

i-Q142, i-Q132

Easy, convenient and reliable blood pressure measurement as never seen before

The OMRON SpotArm™ BPMs are truly comfortable and completely automatic blood pressure monitors. They offer reliable and exact upper arm measurement without having to worry about the cuff. Simply place your arm into the SpotArm™ monitor, switch it on and press start. A few moments later, your blood pressure and pulse values can be viewed in large characters on a 3-fold display. OMRON's innovative technology and highest quality is built-in into every OMRON SpotArm™ Blood Pressure Monitor

The SpotArm BPMs also allow for tracking of Morning Hypertension, a well-known health risk factor. They are fully equipped with OMRON's Intellisense™ technology and irregular heartbeat detection. Several averaging modes plus 84 memories for each of 2 users relieve the user from manual calculations and daily documentation.



HEM-1040-E



HEM-1010-E

- » For arm circumferences from 22–42 cm
- » PC/printer port connector (for optional printer and PC connection)
- » AC powered (Adapter included)

- » For arm circumferences from 17–32 cm
- » Includes batteries and AC adapter

Blood pressure monitors Accessories



- » OMRON SS • 12-18 cm
- » OMRON S • 17-22 cm
- » OMRON M • 22-32 cm
- » OMRON L • 32-42 cm
- » OMRON XL • 42-50 cm

HBP-1300 Cuffs To assure the highest versatility in using the HBP-1300 monitor, OMRON offers these durable and wipeable cuffs in five sizes. It is essential to always use the appropriate size on your patients in order to achieve the highest accuracy.

HBP-1300 Extended Tube

Using the longer replacement tube (either 3.5 m or 1.5m length) allows your patient to sit back and relax.



HEM-907 Stand Take the OMRON HEM-907 easily from patient to patient while storing all accessories in the attached basket using the mobile stand.